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### INSTRUCTIONS FOR USE

#### DESCRIPTION

This device is an operating anoscope for dearterializing hemorrhoidopathy, made of polypropylene. The product is sterile, single use.

#### INDICATIONS

The device is conceived for a temporary use (< 60') in the surgical treatment of hemorrhoidal disease, mainly in 3° degree, by performing a dearterializing hemorrhoidopathy. The device must be exclusively used by experienced physicians with the appropriate surgical techniques.

#### WARNINGS AND PRECAUTIONS

1. Verify that the introducer is properly inserted in the anoscope;
2. Insert the device into the anal canal by placing the handle towards the intergluteal furrow;
3. The guide point on the rotating ferrule indicates the position of the operating window;
4. It is advisable to place the operating window at the biggest hemorrhoidal pedicle;
5. Remove the introducer;
6. Put the stitches on the mucous membrane that protrudes inside the operating window;
7. Insert the introducer;
8. Turn the introducer maneuvering the operating window to the next position;
9. Repeat the previous operations of the points 5 – 6 – 7 – 8 for the following positions;
10. At the end of remove the device.

Use absorbable sutures with a minimum absorption time of 60 days. If a device which has not yet been used comes into contact with instruments or other surgical material used on a patient, or if it has been contaminated with organic fluids, dispose of it in order to prevent the risk of transmission of a viral infection. The surgeon is responsible of applying the standards of good practice in case of the device must be used in presence of pre-existing infections. The device should be stored in a dry place, away from heat sources. The federal law of the United States allows the sale of this product only to physicians or on the order of a physician.

#### STERILIZATION

Device sterilized by Ethylene Oxide (EtO).

The content is sterile only if the packaging has not been opened or damaged; don't use if opened or damaged.

The product is single use and it must not be re-sterilized as:

- it must be used on one patient only;
- not reusable or re-sterilizable, because:
  - the structural integrity of the device could be compromised and the device could be damaged with consequent risk of lesions;
  - there could be a risk of contamination of the device and/or infections in the patient and cross infection, including the transmission of infectious diseases.

#### CONTRAINDICATIONS

Hypersensitivity towards the components.

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